

Introduction

Drugs are most commonly given to the body by the oral route of administration, which leads to greater patient compliance because of the ease of administration. Using modern techniques of drug discovery, chemists can synthesize large numbers of molecules in a short time, but they often create molecules that are poorly absorbed. It is important that molecules with low potential for gastro-intestinal absorption be identified as soon as possible, and discontinued from development.

In vivo measurements of absorption are difficult and slow, and not suitable for screening for absorption potential during drug discovery. A number of physicochemical properties (lipophilicity, solubility, size, hydrogen-bonding characteristics) influence the ability of a molecule to permeate through a membrane by passive diffusion. Moreover, permeability, solubility and lipophilicity are pH-dependent for ionisable drugs, which represent about 70% of the molecules listed in the World Drug Index. Therefore, knowledge of the ionisation constant, K_a , is important. This paper details a method for the rapid measurement of pK_a values.

Materials and Methods

The ProfilerSGA (Figure 1) has been developed to address the need for 'high throughput' measurements of ionisation (pK_a) constants. Large numbers of pK_a measurements can be made very quickly, using small amounts of samples. A patented procedure was devised by creating a linear pH gradient using mixed buffers, injecting the sample, and measuring UV absorbance.



Figure 1. ProfilerSGA (Spectrum Gradient Analyser)

Figure 2 shows a schematic of the instrument. Samples are presented as 10 mM solutions prepared in DMSO in 96-well plates. A flowing pH gradient is produced by mixing two buffer solutions containing mixtures of weak acids and bases that do not absorb significantly in the UV above 250 nm. 20 μ L of each sample is diluted with water and then injected directly into the flowing gradient, which then passes through a diode array spectrophotometer measuring in the UV wavelength range. The buffer has been formulated so that its acid-base titration curve is linear over a wide pH range, such that the pH of the gradient is a linear function of time. The solution pH in the measurement flow cell is therefore proportional to the time lapsed since the start of gradient generation. The sample's K_a values are calculated from the change in UV absorbance at multiple wavelengths as a function of pH.




Figure 2. ProfilerSGA schematic

Initially, acidic buffer (pH 2) is passed through the flow-cell of the diode array spectrophotometer (1 mL/min) and diluted sample is introduced at a constant flow rate (0.2 mL/min). The pH gradient is then changed from pH 2 to pH 12 in a period of 90 minutes and controlled by two computer controlled syringe dispensers that are inversely varied in speed whilst maintaining a constant total flow-rate (1 mL/min). Finally, basic buffer (pH 12) and sample are passed through the flow-cell for a short time before the flow-cell is washed ready for the next sample. Alternatively, the pH gradient is run in the opposite direction from pH 12 to pH 2. UV spectra are recorded every second during this time. The whole procedure of sample dilution, injection and flushing takes about 15 minutes.

The pH gradient was calibrated by injecting standards (at concentrations of 5×10^{-5} - 5×10^{-6} M) with known pK_a values and strong UV absorbance. Two solutions were prepared, the first containing two standards mixed together, the second containing another three standards. Each standard solution was injected into the gradient once before a (96 well) plate of samples and again at the end of the sample plate. The UV absorbance of each standard mixture was measured at 30 wavelengths, and the absorbance values were differentiated with respect to time. The positions of the peak maxima on the time axis indicate the times when the greatest change in UV absorbance at all wavelengths is occurring, which correspond to the pK_a values of the standards. Hence, a linear calibration graph may be plotted of pK_a vs. time of maximum peak height (Figure 3a and 3b).




Figure 3a Derivative Spectra of ProfilerSGA Standards. (Peaks at 50, 58, 77, 99 and 104 secs)





Figure 3b Calibration plot of ProfilerSGA

Peaks for unknown samples can be read directly from this calibration to determine their pK_a values, where samples contain single non-overlapping pK_a values.

In the cases of samples having multiple overlapping pK_a values it is necessary to use a calculation technique known as Target Factor Analysis (TFA). To summarize, spectral data consists of absorbance values measured at different wavelengths and pH values. This can be expressed as a matrix **A**. Invoking Beer's law (**A** = **C****E**), **A** is related to species concentration and absorbivity as follows: the **C** matrix represents the concentration of each independent light-absorbing species of the sample as a function of pH and species and the **E** matrix represents the molar absorptivity of each species (with the inclusion of the optical path length) as a function of pH and wavelength. The **A** matrix based on measured values, but at the start of the calculation, the concentrations and absorptivities in the **C** and **E** matrices are based on guessed values. During the calculation, the **A** matrix is deconvoluted using TFA, in which values of concentration and absorptivity are proposed iteratively until **A** - **C****E** tend to a minimum. When the minimum occurs, the pK_a values required to calculate the distribution of species shown in the **C** matrix are assumed to be the correct pK_a values.

Results and Discussion

Figure 4 shows ProfilerSGA values for 25 compounds plotted against the corresponding literature values measured by pH-metric or traditional UV methods. A linear regression line has been fitted to the data with an R^2 value of 0.996, demonstrating the validity of the new technique. Moreover, the plot contains 41 pK_a values indicating that many of the compounds have two or more ionizable groups.



A=acetylsalicylic acid, B=acyclovir, C=alendazole, D=amiloride, E=clonidine, F=diazepam, G=enrofloxacin, H=famotidine, I=flumequine, J=furosemide, K=hexobarbital, L=homidium, M=hydrochlorothiazide, N=labetalol, O=nicotinic acid, P=niflumic acid, Q=nitrazepam, R=oxitetracycline, S=phenazopyridine, T=phenobarbital, U=quinine, V=sulfacetamide, W=tetracycline, X=uracil, Y=warfarin.

Figure 4. Validation of ProfilerSGA

Warfarin is an example of a compound with a single ionizable group and its pK_a value can be clearly read from the 1st derivative plot of absorbance versus pH (figure 5).




Figure 5. 1st derivative plot of absorbance vs pH for Warfarin

Labetalol, on the other hand, is an example of a compound with two highly overlapping pK_a values. The phenol group causes large change in absorbance and obscures the small absorbance change attributed to the amine group (figure 6a and 6b). In this case, TFA is used to determine both pK_a values (figure 6c and 6d).




Figure 6a. & 6b. Absorbance and 1st derivative plot for labetalol




Figure 6c. & 6d. Distribution of species showing pK_a values at pH 7.4 and 9.4 and molar absorptivity profile. Graph Symbols: HX^+ (x), HX (+), X^- (0).

Conclusions

The ProfilerSGA offers the ability to deploy a 'high throughput' approach to the measurement of pK_a values. The instrument works with small quantities of sample taken from 96 well plates. Analysis times are about 5 minutes per sample thereby providing a throughput of several hundred compounds a day.

References: Hill, A.; Bevan, C.; Reynolds, D. Analytical method and apparatus therefor. International Patent Application WO99/13328, 1999.

Box, K.; Bevan, C.; Comer, J.; Hill, A.; Allen, R.; Reynolds, D. Anal. Chem. 2003, 75(4) 883-892.

Literature pK_a values from various sources: Compounds E, F, J, M, N, O, V: Tam, K. Y.; Takacs-Novak, K. Anal. Chim. Acta 2001, 434, 157-167.

Compounds A, I, P, Q, U: Sirius Technical Application Notes, 1.; Sirius Analytical Instruments Ltd.: Forest Row, RH185DW, UK, 1994.

Compound C: Takacs-Novak, K.; Box, K. J.; Avdeef, A. Int. J. Pharm. 1997, 151(2) 235-248.

Compound X: Sirius Technical Application Notes, 2.; Sirius Analytical Instruments Ltd.: Forest Row, RH185DW, UK, 1995.

Compounds B, D, G: Sammes, P. G.; Drayton, C. J.; Hansch, C.; Taylor, J. B. Comprehensive Medicinal Chemistry; Pergamon Press: Oxford, 1990; Vol. 6.

Compound Y: Avdeef, A.; Box, K. J.; Comer, J. E. A.; Hibbert, C.; Tam, K. Y. Pharm. Res. 1998, 15(2), 209-215.

Compound H: Avdeef, A.; Berger, C. M.; Brownell, C. Pharm. Res. 2000, 17(1), 85-89.

Compounds R, W: Stephens, C. R.; Murai, K.; Brunings, K. J.; Woodward, R. B. J. Am. Chem. Soc. 1956, 78, 4155-4157.

Compound S: Avdeef, A. Curr. Top. Med. Chem. 2001, 1, 277-351.

Compound G: Sirius, unpublished pH-metric results.

Compound L: Sirius, unpublished pH-UV result.