



Certificate of Conformity

Product Name: Acidic Buffer Reagent

Part Number: 1404001

Batch Number: 512457

Method of Analysis: pH-metric & spectrophotometric titration


DataSet Numbers: 03D-0207, 03D-0208, 03D-0209, 03D-0210

Linearity Plot for pH3 - 11:

(Specification = r^2 to be ≥ 0.9900)

Average Initial pH = 2.7 ± 0.0

r^2 range = 0.9987 – 0.9989




UV Profile:

(Specification = Absorbance to be ≤ 0.05 between 250nm and 500nm)

DataSet 03D-0208

Wavelength	Absorbance
250nm	0.01
260nm	0.01
280nm	0.01
300nm	0.01
350nm	0.00
400nm	0.00
450nm	0.00
500nm	0.00



This product has been tested by the Quality Control Laboratory of Sirius Analytical Instruments Limited and conforms to the specifications described above.

Approved:

.....
Carolyn Hibbert, Quality Assurance Co-ordinator

Certificate Date:09 April 2003.....

