

Product Name: SGA Calibration Standard 3 **Part Number:** 1404006


Batch Number: C3SGA003

Method of Analysis: Spectrophotometric titration
DataSet Numbers: 04B-1204

pK_a values:

Equation	Measured value	Specification value
$X + H = XH$	9.2	9.2 ± 0.2
$XH + H = XH_2$	5.0	5.0 ± 0.2

UV profile:



This product has been tested by the Quality Control Laboratory of Sirius Analytical Instruments Limited and conforms to the specifications described above.

Approved:..... *Carolyn Hibbert*
 Carolyn Hibbert, Quality Assurance Co-ordinator

Certificate Date:..... 12 February 2004.....

C3SGA003.doc